


THE EUROPEAN MALIGNANT HYPERTHERMIA GROUP

Ryanodine test - optional and standardized for the EMHG

<p>Ryanodine</p> <p>Stock solution</p> <p>Storage</p> <p>Test concentration</p> <p>General guidelines</p>	<p>High purity (Calbiochem).</p> <p>100 micromol litre⁻¹ in distilled water.</p> <p>Refrigerator at 4°C.</p> <p>Final bath concentration of 1 micromol litre⁻¹.</p> <p>The ryanodine test is performed similarly to other contracture tests described in the EMHG protocol:</p> <ol style="list-style-type: none"> 1. A fresh muscle specimen must be used. 2. The specimen is electrically stimulated with a 1 ms supra-maximal stimulus at a frequency of 0.2 Hz. 3. The test is performed at optimal length (l_o). 4. Baseline must not vary more than 2 mN (0.2 g) within a 10 min period before the addition of ryanodine. 5. A single bolus dose technique is used to reach a final bath concentration of 1 micromol litre⁻¹. If the tissue bath is continuously perfused with Krebs solution, either the perfusion is stopped or a Krebs-ryanodine solution is perfused. The final bath concentration of ryanodine must be reached within 1 min.
<p>Test parameters</p>	<ol style="list-style-type: none"> 1. Similar to other tests: weight (mg) and length (mm) of specimen, twitch height (mN), maximum preload (mN), predrug baseline (mN), lowpoint baseline (mN). 2. Specific parameters: onset time (OT) of contracture (min), time for development of a 10 mN (1 g) contracture (10T) (min), contracture velocity (dT = 10T-OT) (min), and optionally: time for development of contracture maximum (T_{max}) (min), maximum contracture (mN).
<p>Quality control</p>	<p>The concentration of ryanodine in the tissue bath must be periodically checked by HPLC.</p>



In this example, predrug baseline = low point baseline.